


CQ-01717-D2O-DEPT135-271009-298 K-500 MHz
 DEPT135-mmcd-clem D2O /opt/topspin ravi {1 D2 - 138}

54.26
 53.26
 53.16
 33.21
 28.24
 24.96
 24.86
 24.59
 18.99

— -2.82



NAME CQ-01717-D2O-271009
 EXPNO 5
 PROCNO 1
 Date_ 20091101
 Time 2.37
 INSTRUM spect
 PROBHD 5 mm CPTXO 13C
 PULPROG dept135
 TD 65536
 SOLVENT D2O
 NS 256
 DS 16
 SWH 29761.904 Hz
 FIDRES 0.454131 Hz
 AQ 1.1010548 sec
 RG 203
 DW 16.800 usec
 DE 30.00 usec
 TE 298.0 K
 CNST2 145.0000000
 D1 2.00000000 sec
 D2 0.00344828 sec
 D12 0.00002000 sec
 TD0 1

===== CHANNEL f1 =====
 NUC1 13C
 P1 10.00 usec
 P2 20.00 usec
 PL1 3.00 dB
 PL1W 36.25705338 W
 SFO1 125.6974365 MHz

===== CHANNEL f2 =====
 CPDPRG2 waltz16
 NUC2 1H
 P3 15.70 usec
 P4 31.40 usec
 PCPD2 100.00 usec
 PL2 1.25 dB
 PL12 17.33 dB
 PL2W 19.96668053 W
 PL12W 0.49238619 W
 SFO2 499.8419994 MHz
 SI 32768
 SF 125.6848680 MHz
 WDW EM
 SSB 0
 LB 1.00 Hz
 GB 0
 PC 1.40

